

Q1. How can I uninstall the Mirascan driver from Windows XP

Solution: If you cannot uninstall the Mirascan driver in Windows XP, please follow these steps:

Step 1: Load the Windows Registry Editor by running REGEDIT

Step 2: Press F3 and search for MIRASCAN

Step 3: For each matching value or folder that contains MIRASCAN, press delete and remove the entry. **DO NOT** delete entries that you are unsure of, or that do not contain MIRASCAN

Step 4: Press F3 again to search for the next value and continue to delete all MIRASCAN entries until the search finds no matching values

Step 5: Restart Windows XP and reinstall the driver from the www.benq.com website